

Solution Requirements

Diabetic Retinopathy Detection System

FUNCTIONAL REQUIREMENTS

1. User Management

FR-1.1: User Registration - Allow new users to create accounts with name, email, password - Validate email format and password match - Verify email uniqueness in database - Priority: MUST HAVE

FR-1.2: User Login - Authenticate users with email and password - Create secure session upon successful login - Priority: MUST HAVE

FR-1.3: User Logout - End user session securely - Clear session data - Priority: MUST HAVE

2. Image Upload & Processing

FR-2.1: Image Upload - Accept retinal fundus images (PNG, JPG, JPEG) - Validate file format and size (max 16MB) - Secure filename handling - Priority: MUST HAVE

FR-2.2: Image Preprocessing - Resize images to 299x299 pixels - Normalize pixel values (0-1 range) - Convert to model input format - Priority: MUST HAVE

3. DR Classification

FR-3.1: DR Prediction - Classify DR severity into 5 classes: 1. No_DR, 2. Mild, 3. Moderate, 4. Severe, 5. Proliferate_DR - Load trained Xception model - Perform inference and extract probabilities - Priority: MUST HAVE

FR-3.2: Confidence Score Display - Display prediction confidence as percentage (0-100%) - Show probabilities for all 5 classes - Priority: MUST HAVE

4. Results Management

FR-4.1: Prediction Display - Show DR classification, confidence, uploaded image, timestamp - User-friendly result presentation - Priority: MUST HAVE

FR-4.2: Prediction History Storage - Store predictions in database with user, result, confidence, timestamp - Enable future retrieval - Priority: SHOULD HAVE

5. Navigation & Interface

FR-5.1: Responsive Design - Adapt interface to desktop, tablet, mobile - Touch-friendly on mobile devices - Priority: SHOULD HAVE

NON-FUNCTIONAL REQUIREMENTS

Performance

NFR-1.1: Prediction Speed - ≤ 5 seconds per image (HIGH)

NFR-1.2: Page Load Time - < 3 seconds (MEDIUM)

NFR-1.3: Concurrent Users - Support minimum 10 users (MEDIUM)

Accuracy

NFR-2.1: Classification Accuracy - $\geq 85\%$ overall (CRITICAL)

NFR-2.2: Class-wise Performance - F1-score ≥ 0.80 per class (HIGH)

NFR-2.3: False Negative Rate - Minimize for severe cases (CRITICAL)

Reliability

NFR-3.1: System Availability - 99% uptime (HIGH)

NFR-3.2: Error Handling - Graceful with user-friendly messages (HIGH)

NFR-3.3: Data Integrity - No data loss, regular backups (HIGH)

Security

NFR-4.1: Authentication - Secure session-based auth (CRITICAL)

NFR-4.2: Data Privacy - Protected user data, HIPAA considerations (CRITICAL)

NFR-4.3: Secure File Upload - Validate and sanitize files (HIGH)

NFR-4.4: Password Security - Bcrypt hashing for production (CRITICAL)

Usability

NFR-5.1: User Interface - Intuitive for non-technical users (HIGH)

NFR-5.2: Accessibility - WCAG 2.1 guidelines (MEDIUM)

NFR-5.3: Error Messages - Clear, actionable, non-technical (MEDIUM)

Compatibility

NFR-6.1: Browser Support - Chrome, Firefox, Safari, Edge (latest 2 versions) (HIGH)

NFR-6.2: Operating System - Windows, macOS, Linux, iOS, Android (HIGH)

NFR-6.3: Screen Resolutions - 320px to 1920px+ (MEDIUM)

DATA REQUIREMENTS

Training Data

- Minimum 3,000 labeled fundus images (Actual: 3,662)
- All 5 DR severity levels represented
- High-quality PNG images

User Data

- User information: Name, email, password, registration date
 - Prediction records: User, prediction, confidence, timestamp, image name
 - Storage: JSON documents in IBM Cloudant
-

INTERFACE REQUIREMENTS

User Interface

- Registration Form: Name, Email, Password, Confirm Password
- Login Form: Email, Password
- Prediction Interface: File upload, submit button, result display

API Interface

- Model Inference API: Input (image array) → Output (class probabilities)
 - Database API: CRUD operations via Cloudant REST API
-

CONSTRAINTS

Technical: Python 3.8+, TensorFlow 2.x, JavaScript-enabled browser, internet connection

Resource: 4-6 week timeline, 1-3 developers, minimal budget

Regulatory: Medical device regulations, data privacy laws (HIPAA, GDPR)

ASSUMPTIONS

1. Users have basic computer literacy
 2. Retinal images are properly captured
 3. Internet connectivity available
 4. Modern web browser installed
 5. IBM Cloud account for Cloudant
 6. Training data properly labeled
-

ACCEPTANCE CRITERIA

- ✓ Users can register and login
- ✓ Users can upload retinal images
- ✓ System classifies DR with $\geq 85\%$ accuracy
- ✓ Predictions display within 5 seconds
- ✓ Results show class and confidence
- ✓ Predictions stored in database
- ✓ Interface is responsive

- ✓ System handles errors gracefully
 - ✓ Security measures implemented
 - ✓ Documentation complete
-

OUT OF SCOPE (Future Versions)

Mobile app, batch processing, PDF reports, email notifications, EHR integration, multi-language support, heatmap visualization, admin dashboard, analytics